



6-BENZYLAMINOPURINE RIBOSIDE

CAS Number: 4294-16-0

CERTIFICATE OF ANALYSIS

Molecular Weight: 357.37

A white, microcrystalline solid.

M.P. 170.0 – 171.0 °C.


UV: water, λ_{\max} 268.5 nm

HPLC purity, 268 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches $[M]^+ = 358.16$ amu.

^1H NMR (DMSO): in accordance

Synonym: BAR/BAPR, N⁶-BENZYLADENOSINE



BioApex s.r.o.

Technology Park of Palacký University

Šlechtitelů 21

783 71 Olomouc - Holic

Czech Republic

Phone: +420 585 631 420

Cell phone: +420 737 815 152

E-mail: bioapex@atlas.cz