

6-BENZYLAMINOPURINE

CAS Number: 1214-39-7

CERTIFICATE OF ANALYSIS

Molecular Weight: 225.25

A white, microcrystalline solid.

M.P. 230.0 – 233.0 °C.


UV: water, λ_{\max} 269.3 nm

HPLC purity, 269 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches $[M]^+ = 226.06$ amu.

^1H NMR (DMSO): in accordance

Synonym: BA/BAP, N⁶-BENZYLADENINE



BioApex s.r.o.

Technology Park of Palacký University

Šlechtitelů 21

783 71 Olomouc - Holice

Czech Republic

Phone: +420 585 631 420

Cell phone: +420 737 815 152

E-mail: bioapex@atlas.cz