

DIHYDROZEATIN
CAS Number: 14894-18-9

CERTIFICATE OF ANALYSIS

Molecular Weight: 221.26

A colourless, amorphous solid.

M.P.: 169 – 170 °C


UV: water, λ_{\max} 269.3 nm

HPLC purity, 268 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches $[M]^+ = 222.18$ amu.

^1H NMR (DMSO): in accordance

Synonym: DZ/DHZ, 6-(4-hydroxy-3-methylbutylamino)purine



BioApex s.r.o.
Technology Park of Palacký University
Šlechtitelů 21
783 71 Olomouc - Holic
Czech Republic
Phone: +420 585 631 420
Cell phone: +420 737 815 152
E-mail: bioapex@atlas.cz