

N⁶-(Δ²-ISOPENTENYL)ADENINE

CAS Number: 2365-40-4

CERTIFICATE OF ANALYSIS

Molecular Weight: 203.24

Colourless, microcrystalline solid.

M.P. 210.0 – 211.3 °C.


UV: water, λ_{max} 269 nm

HPLC purity, 268 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches [M]⁺ = 204.2.

¹H NMR (DMSO): in accordance

Synonym: 2iP/iP, 6-(γ,γ-Dimethylallylamino)purine



BioApex s.r.o.

Technology Park of Palacký University

Šlechtitelů 21

783 71 Olomouc - Holice

Czech Republic

Phone: +420 585 631 420

Cell phone: +420 737 815 152

E-mail: bioapex@atlas.cz