

N⁶-(2-ISOPENTENYL)ADENOSINE

CAS Number: 7724-76-7

CERTIFICATE OF ANALYSIS

Molecular Weight: 335.37

Colourless, crystalline solid.

M.P. 148.0 – 150.0 °C.


UV: water, λ_{\max} 269.3 nm

HPLC purity, 268 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches $[M]^+ = 336.2$.

¹H NMR (DMSO): in accordance

Synonym: iPR, 6-(γ,γ -dimethylallylamino)purine riboside



BioApex s.r.o.

Technology Park of Palacký University

Šlechtitelů 21

783 71 Olomouc - Holice

Czech Republic

Phone: +420 585 631 420

Cell phone: +420 737 815 152

E-mail: bioapex@atlas.cz