



MYOSEVERIN

CAS Number: 267402-71-1

CERTIFICATE OF ANALYSIS

Molecular Weight: 432.10

A colourless solid.

M.P.: 109-110 °C


UV: 90% EtOH /0.1 M HCl; λ_{\max} 258(sh) , 284, 296 nm

HPLC purity, 235.0 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches $[M]^+ = 433.2$ amu.

^1H NMR (DMSO): in accordance

Synonym: 9-isopropyl-2,6-bis(4-methoxybenzylamino)purine



BioApex s.r.o.

Technology Park of Palacký University

Šlechtitelů 21

783 71 Olomouc - Holice

Czech Republic

Phone: +420 585 631 420

Cell phone: +420 737 815 152

E-mail: bioapex@atlas.cz