

R-ROSCOVITINE
CAS Number: 186692-46-6

CERTIFICATE OF ANALYSIS

Molecular Weight: 354. 45

A colourless solid.

M.P.: 136-139°C


UV: 90% EtOH /0.1 M HCl; λ_{\max} 294 nm

HPLC purity, 235.0 nm, 99+%.

MASS SPECTROMETRY: Passes the ESI(+) test; matches $[M]^+ = 355.3$ amu.

^1H NMR (DMSO): in accordance

Synonym: ROSC, 6-benzylamino-[2-(1(R)-(hydroxymethyl)propyl)amino-9-isopropylpurine



BioApex s.r.o.
Technology Park of Palacký University
Šlechtitelů 21
783 71 Olomouc - Holice
Czech Republic
Phone: +420 585 631 420
Cell phone: +420 737 815 152
E-mail: bioapex@atlas.cz